



Release Notes: F-950

Version 1.9.5.4

Release Date: August 9, 2017

Additions

- Support for screen capture using G.A.S software

Changes

- External KMNO4 chamber required during a “Set Zero” procedure (does not use internal chamber)
- Set zero now takes 600 seconds (previously it was 300)
- Bug fixes

Version 1.9.4.1

Release Date: March 12, 2017

Additions

- Passcode to enter calibration menu enabled (SETUP>CALIBRATION) as well as the Set Zero menu (SETUP>SET ZERO).

Changes

- Menus related to calibration got rearranged
-

Version 1.9.3.3

Release Date: March 7, 2017

Initial release for Hardware V4 of the F-950
